

Gellért Sipos, László Lengyel, György Dormán, László Kocsis, Ferenc Darvas and Richard Jones

ThalesNano Inc., Záhony u. 7., H-1031, Budapest, Hungary


richard.jones@thalesnano.com, www.thalesnano.com

Introduction

Gould-Jacobs type intramolecular thermal cyclisations were previously reported about in a continuous flow reactor at high temperatures (300-360 °C) and pressure (100-160 bar) in liquid phase [1]. The regioselectivity of the ring closure is dependent on the nature and position of the substituents often leading to a mixture of products. We investigated the regiochemical outcome of the cyclisation of substituted amino-pyridine substructures under various conditions in liquid as well as in gas phase, by applying a flow pyrolysis apparatus under high vacuum in solvent-free conditions. Flash vacuum pyrolysis (FVP) allows a rapid exposure to high temperature (200 - 900 °C), which often favours the formation of one regioisomer. Here we compare the outcome of the thermal cyclisation of various unsaturated diesters or analogous substructures (derived from condensation with ethoxymethylenemalonate) performed in both systems.

Thermal cyclisation – the Gould-Jacobs reaction

Alkylidene β-diester or analogous substructures derived from condensation with Meldrum's acid or electron-withdrawing group containing esters (e.g. malonic ester or cyanoacetic acid esters) readily undergo pericyclic annulation reactions at high temperatures, if they are located in a favourable position relative to aromatic or unsaturated systems (Scheme 1.).



Scheme 1.: The Gould-Jacobs reaction

Typical procedures

The thermal cyclisation of the above-mentioned types is normally carried out in high boiling point solvents (e.g. diphenyl ether). In several cases, the long exposure to high temperature induces side-reactions. The work-up can therefore be difficult. In high temperature/pressure continuous flow systems, such as the Phoenix Flow Reactor (Scheme 2.), low boiling point solvents can be used with high flow rates to avoid side reactions or decomposition. The use of microwave or gas phase techniques can be beneficial as well.

Solution phase flow pyrolysis (SPFP) in the Phoenix Flow Reactor



Scheme 2.: Phoenix Flow Reactor

Flash vacuum pyrolysis (FVP)

Pyrolysis is a thermochemical treatment of organic materials at elevated temperatures in gas phase (~200-900°C). FVP is a technique which involves a vacuum distillation or sublimation of a substrate through an empty or filled hot tube with the products afterwards collected in a cold trap. The major elements of the device are: Evaporation unit, Reactor unit, Condenser unit, Vacuum unit. Typical residence time: ms – s.

Regioselectivity of the Gould-Jacobs reaction

The regioselectivity of the Gould-Jacobs type reactions is another important issue, which is normally governed by steric and electronic factors and often leads to a mixture of products. For example, 2-aminopyridine-derived alkylidene esters may be applied for the synthesis of various pyrido[1,2-a]pyrimidin-4-ones and 1,8-naphthridinones depending on the substituents on the pyridine ring and the other substituents [2].

Herein, we present the study of the ring closure of *N*-(*x*-methyl-2-pyridyl)aminomethylenemalonate derivatives (I.a-c, VI.) using SPFP and FVP method, as well.

Results and discussion

The gas phase pyrolysis of I.a-c at 450°C exclusively led to the formation of pyridopyrimidones (III.a-c), without other side products or starting material observed during the course of the reactions. However, applying the solution phase flow pyrolysis method (for the above mentioned substrates) resulted in product mixtures. In the case of I.a the main product was III.a with 60% conversion, the decarboxylated V.a was formed with 21%. Furthermore, 19% of unidentified products also came about. It should be noted that regarding the ring closure of I.a, the formation of naphthridone is not possible due to the position of the methyl group. When the reaction was attempted on substrate I.b, two main products: III.b. and V.b were formed approximately in a 6:4 ratio respectively. Finally, III.c and V.c developed in a nearly 1:1 ratio during the solution phase flow pyrolysis of I.c. Moreover, in the latter case we observed the formation of II.c (1%). Interesting to note that if we carried out the cyclisation reactions of I.a-c using higher flow rates (therefore, lower residence time), we observed higher amounts of the decarboxylated compounds (V.a-c). Since at the moment we cannot give a logical reason for these results, we would like to extend our studies to understand this phenomenon.

Summary

We examined the ring closure reaction of all the *N*-(*x*-methyl-2-pyridyl)aminomethylenemalonate derivatives in solution and in gas phase as well. When the methyl group was in the 3-, 4- or 5-positions, the results of the solution and gas phase methods were roughly the same. The only difference being the rate of decarboxylation in solution phase depending on the applied flow rate.

However, important differences were observed during the cyclisation of *N*-(6-methyl-2-pyridyl)aminomethylenemalonate (VI.): While the FVP reaction led exclusively to the pyridopyrimidone, during the SPFP, a complex product mixture was formed. The expected naphthridone structures VIII. and IX. contribute 18% and 14% respectively. Notably, the main product, X. made 49%. The latter can be converted to IX. by applying further heating on the substrate.





Table 1.: Thermal cyclisation of I.a-c

Method	Substrate	Conversions (%)				
		II.	III.	IV.	V.	Other
FVP	I.a	-	100	-	0	0
FVP	I.b	0	100	0	0	0
FVP	I.c	0	100	0	0	0
SPFP	I.a	-	60	-	21	19
SPFP	I.b	0	57	0	38	5
SPFP	I.c	1	47	0	46	6

Conditions for SPFP: 350 °C, 100 bar, 4 ml/min, 4 ml loop volume
Conditions for FVP: 450 °C, 10⁻³ mbar, 30 cm tube length

Pyrolysis of *N*-(6-methyl-2-pyridyl)aminomethylenemalonate (VI.)

According to literature data, the thermal ring closure of VI. in batch (with microwave or conventional heating) exclusively leads to naphthridone (VIII.) due to the steric hindrance of the methyl group [3]. However, the pyridopyrimidone (VII.) can be synthesized in an acid catalysed reaction [4].



Scheme 4.: Pyrolysis of VI.

We found that during the gas phase thermal cyclisation of VI. the pyridopyrimidone (VII.) was formed as an only product. Unexpectedly, when we carried out the reaction in the X-Cube Flash, compound X. was the main product (49%). We also observed the formation of VIII. and IX., but we were not able to detect VII. in the reaction mixture.

The above-mentioned formation of X. means that decarboxylation occurs as the first reaction step. After that, via a ring closure, it may turn into IX.. In order to confirm this idea, we pumped through the previously isolated X. under the following conditions: 350°C, 100 bar, 4 ml/min, 4 ml loop. The analytical results proved X.'s rearrangement into IX.

Table 2.: Thermal cyclisation of VI.

Method	T / °C	Conversions (%)					
		VI.	VII.	VIII.	IX.	X.	Other
FVP	250	100	0	0	0	0	0
FVP	350	60	40	0	0	0	0
FVP	400	0	100	0	0	0	0
FVP	500	0	100	0	0	0	0
SPFP	350	7	0	18	14	49	12

Experimental

Typical SPFP reaction: a 9 mM or a 0.1 M CH₃CN solution of I.a was pumped through the X-Cube Flash™, under the following conditions: 350 °C, 100 bar, 4 ml/min, 4 ml loop volume. The product was analysed with LC-MS and NMR. (If it was necessary column chromatography was also carried out.)

Typical FVP reaction: 100 mg of I.a was pyrolysed under the following conditions: 30 cm tube length, T_{urnace} = 450°C, T_{inlet} = 130-140°C, p = 10⁻³ mbar, pyrolysis time = 1 h. The product was washed out from the collection tube with DCM, then the solvent was evaporated. Yield: 98%. The crude product was analysed with LC-MS and NMR.