

NEW
Cooled Carousel
6 Place Now Available



Introductory Leaflet

Carousel Work-Up Station™



Rapid parallel work-up
using standard 70ml
columns for efficient
filtration, phase
separation or SPE...




Carousel Work-Up Station

Rapid parallel work-up using standard 70ml columns for efficient filtration, phase separation or SPE...

- Compatible with 12 x Carousel Reaction Tubes or standard 1 inch boiling tubes
- Use under gravity or with new SpeediFlow flow booster
- Accepts 70ml columns; including 20um Filtration, Phase Separation, Silica, SCX, Aminopropyl and SAX Acetate
- Simple and convenient the Carousel Work-Up Station will reduce post synthesis bottlenecks...
- Identical aluminium racks are easy to load with tubes or columns
- Racks are designed to lock together for ease of transportation
- Optional 3rd rack can be stacked above the 2nd rack for sequential purification




Personal synthesis stations for parallel solution phase chemistry and solid supported reagent based synthesis



Cost effective low temperature parallel synthesis from ambient to -70°C. Accepts 6 x 100ml, 170ml or 250ml reactors

Cat No	Description	Pk Qty
RR99830	Carousel Work-Up Station (Two racks and SpeediFlow)	1
Carousel Work-Up Station Components & Accessories		
RR99831	Carousel Work-Up Rack	1
RR99850	SpeediFlow	1
RR99860	Replacement Viton O-ring - 1 in OD	5
RR99861	Replacement Perlast O-ring - 1 in OD	2
RR99870	Silica 10g/70ml	16
RR99871	Aminopropyl 10g/70ml	16
RR99872	SCX 10g/70ml	16
RR99873	SAX Acetate 10g/70ml	16
RR99874	20um PE Filtration 70ml	50
RR99875	Phase Separation 70ml	50
RR99840	SpeediFlow Gas Connection Tubing - 5m length	1



Save time and money by buying the popular Carousel Work-Up 'Systems', which include all the important components and accessories required for efficient operation.

Cat No	Description	Pk Qty
Carousel Work-Up Column Pack		
RR99876	Column Pack Comprising of the following: RR99870 - Silica 10g/70ml, pack of 16 RR99871 - Aminopropyl 10g/70ml, pack of 16 RR99872 - SCX 10g/70ml, pack of 16 RR99873 - SAX Acetate 10g/70ml, pack of 16 RR99874 - 20um PE Filtration 70ml, pack of 50 RR99875 - Phase Separation 70ml, pack of 50	1
Carousel Work-Up System 25		
RR99880	System 25 - Carousel Work Up Station Comprising of the following: RR99830 - Carousel Work-Up Station (Two racks and SpeediFlow) RR99874 - 20um PE Filtration 70ml, pack of 50 RR99875 - Phase Separation 70ml, pack of 50 RR98062 - Threaded Glass Reaction Tube 24mm x 150mm, pack of 12 RR99840 - SpeediFlow Gas Connection Tubing - 5m length	1
Carousel Work-Up System 26		
RR99881	System 26 - Carousel Work Up Station Comprising of the following: RR99830 - Carousel Work-Up Station (Two racks and SpeediFlow) RR99831 - Carousel Work-Up Rack RR99870 - Silica 10g/70ml, pack of 16 RR99874 - 20um PE Filtration 70ml, pack of 50 RR99875 - Phase Separation 70ml, pack of 50 RR98062 - Threaded Glass Reaction Tube 24mm x 150mm, pack of 12 RR99840 - SpeediFlow Gas Connection Tubing - 5m length	1

Radleys Discovery Technologies are specialists in parallel chemistry...

As a dynamic organisation Radleys Discovery Technologies are best able to react to the needs of this constantly developing market sector. Radleys Discovery Technologies specific areas of expertise is focused on apparatus and consumables for parallel synthesis, purification and work-up. The essential products that are the cornerstone of your parallel chemistry program.

© Radleys Discovery Technologies 2003