

Frozen in Time Ltd.

Manufacturers of Freeze Drying Machines
and Vacuum Cold traps

Lablyo Freezedrier

UK Manufactured

Bench top freeze drier with new Touch Screen controls. Data logging and display of temperature and vacuum on real time graph. Data downloadable via USB. 20cm diameter x 35cm deep ice condenser chamber with 5kg ice capacity. 3kg ice collection in 24 hours. Condenser temperature -55°C . Mobile temperature probe for product sampling.


Part No's

Lablyo	FIT/LYO/49/1
Vacuum pump for Lablyo	FIT/LYO/59/0

Can be used with the following freeze drying accessories:

Part No's

Acrylic drying chamber	FIT/LYO/01/0
6 height shelf rack	FIT/LYO/03/0
6 heated shelf rack	FIT/LYO/19/0
4 shelf vial closing device	FIT/LYO/22/0
8 port flask manifold	FIT/LYO/13/0
48 port ampoule manifold	FIT/LYO/25/0
Spin Freezer	FIT/LYO/28/0



FROZEN IN TIME Ltd

York Road Industrial Park, Sheriff Hutton, York, YO60 6RZ England.
TEL. 01347 878158. FAX. 01347 878303. E-MAIL. info@freezedriers.com WEBSITE. www.freezedriers.com